

Software Architecture

This document explains the overall design of the diode-characterization system: what it does, how a measurement is configured, how it is run, what each settings file controls, and why there are several "engine" files and how each one is selected and driven.

1. The big idea

This is an **automated lab-bench measurement system** for characterizing diodes. It drives three Keithley instruments over GPIB (via PyVISA) and records the results to CSV, then rolls those CSVs up into a single pass/fail report in the test station's text format.

The three instruments:


Instrument	Role	Driver file
Keithley 4200-SCS	Source-Measure Unit (SMU) — sources V or I, measures V & I	keithley4200.py
Keithley 590	LCR meter — measures capacitance & conductance	keithley590_LCR.py
Keithley 707A	Switching matrix — routes signals to the device pin under test	keithley707a.py

The **central loop** is always the same shape, regardless of measurement type:

```
for each row in the Excel parameter matrix:
    1. set the switching matrix to connect the right pin    (707A)
    2. apply source setpoints and let them settle          (4200 or 590)
    3. measure                                              (4200 or 590)
    4. append a result row
finally:
    safe shutdown (power off sources, open the matrix)
```

What *changes* between measurement types is **what gets sourced and measured in steps 2–3**. That variation is captured by the different engine classes (see §6).

Data flow at a glance



```
results/<name>_<ts>.csv

run_all.py aggregates 4 CSVs + settings/run_info.json → ReportWriter →
results/result_*.txt
```

2. How a measurement is configured

A measurement is configured by **two inputs**:

- 1. **A JSON config** in `settings/` — *how* to measure: instrument addresses, source setpoints, compliance, speed, etc. One JSON file == one measurement type.
- 2. **An Excel workbook** (`settings/ParameterMatrix_DF.xlsx`) — *where* to measure: one row per device pin, giving the switching-matrix routing and the pin label. Each JSON config names which **sheet** of the workbook it reads.

So the JSON answers "what stimulus do I apply and how do I read it back," and the Excel sheet answers "which pins do I iterate over and how do I wire them up." The engine takes the cross product: for every row in the sheet, it runs the stimulus defined in the JSON.

3. The JSON settings files

All configs live in `settings/`. There are four measurement configs plus one run-metadata file.

File	Measurement	Engine selected
<code>guard_leakage.json</code>	Guard ring leakage current	<code>SweepRunner</code>
<code>dark_current.json</code>	Reverse/forward dark current	<code>SweepRunner</code>
<code>series_resistance.json</code>	Series resistance Rs	<code>SeriesResistanceRunner</code>
<code>capacitance.json</code>	Junction capacitance	<code>LCRRunner</code>
<code>run_info.json</code>	Wafer/lot/operator metadata for the report (not a measurement)	—

3.1 Common top-level keys

```
{
  "test_name": "Dark Current",           // used for CSV filename + report
  "excel_sheet": "DarkCurrent",         // which workbook sheet to iterate
  "instruments": {
    "switch_matrix": "GPIB0::18::INSTR", // 707A VISA address
    "smu_mainframe": "GPIB0::17::INSTR", // 4200 VISA address (SMU configs)
    "lcr_meter":     "GPIB0::15::INSTR",  // 590 VISA address (LCR config only)
    "matrix_settling_s": 0.3,             // relay settle delay after routing
    "workbook": "settings/ParameterMatrix_DF.xlsx" // optional; overrides default
  }
}
```

```
// ... then either "smus" or "lcr" (see below)
}
```

3.2 SMU configs (smus block)

Used by `guard_leakage.json`, `dark_current.json`, `series_resistance.json`.

Each entry under "smus" is one SMU channel. **The key encodes the role and the SMU id: "cathode_SMU1"** → role `cathode`, channel `SMU1`. The *role* becomes the CSV column prefix (`Cathode Current A`, etc.).

```
"cathode_SMU1": {
  "sweep_mode": "Voltage in V", // or "Current in A" (default: Voltage in V)
  "compliance_A": 0.1,         // current/voltage compliance limit
  "speed": "Slow",             // 4200 integration speed
  "range": "Auto",             // measurement range
  "average": 1,                 // readings averaged per point
  "sweep_profile": [           // list of (setpoint, hold) steps
    {"voltage": -1.25, "hold_s": 1.0},
    {"voltage": 2.5, "hold_s": 3.0}
  ]
}
```

Each `sweep_profile` is a loop axis — the single most important rule:

`SweepRunner` treats every channel's `sweep_profile` as one axis of a **nested sweep** (Cartesian product). For each Excel row it sweeps every combination of every channel's steps:

- **JSON declaration order = loop nesting order.** The **first** channel in the JSON is the **outermost** loop (varies slowest); the **last** channel is the **innermost** (varies fastest). All channels are applied in declaration order at every combination, each honoring its own `hold_s` — there is no implicit skipping, so the timing is exactly what the JSON declares.
- A channel with a **single** step is simply held at that value across the whole nest, so it contributes one point to its axis. A config where only one channel is multi-step therefore collapses to a plain 1-D sweep of that channel — which is exactly how `guard_leakage.json` and `dark_current.json` behave.
- **Any number** of multi-step channels is allowed; the total row count per pin is the product of the per-channel step counts. `load_measurement_spec()` no longer restricts this. (The `ChannelSpec.is_primary` helper — "more than one step" — still exists and is used by the series-resistance engine to locate its swept channel.)

`sweep_mode` lets a channel source current instead of voltage (used by series resistance, which forces two current setpoints). Despite the field being named `voltage`, it carries the current value in amps when `sweep_mode` is "Current in A".

3.3 LCR config (lcr block)

Used by `capacitance.json`. No SMUs; instead an "lcr" block configures the 590:

```

"lcr": {
  "frequency": 100000000.0,      // measurement frequency (Hz)
  "integration": "10 /s",        // integration / reading rate
  "range": "Auto",
  "trigger": "One-shot, talk",
  "bias_voltage": 0.0            // DC bias applied during the C/G read
}

```

The presence of the "lcr" key is itself the dispatch signal (see §5).

3.4 run_info.json

Pure metadata — wafer id, lot id, operator, device position, temperature, test station. It is **not** a measurement config; it is consumed only by the report writer to fill in the report header and to build the output filename.

4. The Excel parameter matrix

`parameter_matrix.py` reads the workbook. `load_parameter_rows()`:

- Opens the workbook, selects the sheet named by `excel_sheet`.
- **Normalizes column names** (case- and whitespace-insensitive; e.g. `matrixconfig`, `Matrix Config`, `MATRIX CONFIG` all map to `Matrix Config`).
- Requires two columns: **Matrix Config** (the 707A crosspoint routing string) and **Measured Pin** (the pin label). Rows missing either are skipped with a console note.
- Returns a `List[Dict[str, str]]` — one dict per pin.

Each measurement type uses a **different sheet** of the same workbook, so the same physical file holds the pin list for guard leakage, dark current, capacitance, and series resistance independently.

The matrix routing string itself is canonicalized by `normalize_matrix_config()` in `switching_matrix.py` (uppercased, comma/semicolon-tolerant, joined with `;`).

5. How to run

Single measurement — `sweep.py`

```

python sweep.py settings/dark_current.json
python sweep.py settings/guard_leakage.json --dry-run
python sweep.py settings/capacitance.json --limit-rows 3
python sweep.py settings/series_resistance.json --output results/rs.csv

```

Flags:

Flag	Effect
<code>--dry-run</code>	Simulate without touching hardware or sleeping (uses <code>DryRunResourceManager</code> / <code>DryRunPort</code> , which print every driver command).
<code>--limit-rows N</code>	Process only the first N rows of the Excel sheet.
<code>--output PATH</code>	Override the CSV output path (default: <code>results/<TestName>_<timestamp>.csv</code>).

One device, no prober — `run_all.py`

```
python run_all.py
python run_all.py --dry-run --limit-rows 2
```

`run_all.py` runs the four configs **in a fixed order** (guard leakage → dark current → capacitance → series resistance), waiting a few seconds between tests, then feeds the four resulting CSVs plus `run_info.json` into `ReportWriter`. The reusable core is `run_device(run_info, *, dry_run, limit_rows, output_dir)`, which writes all four CSVs and the report into one **per-device folder** `results/wafer<id>/dev<NNN>_x<x>_y<y>/` (fixed, collision-free names). This is the **no-prober** path — `run_all.py` never imports the prober.

Full wafer — `run_wafer.py`

```
python run_wafer.py GPIB0::1::INSTR accretech/examples/DF.mdf
python run_wafer.py GPIB0::1::INSTR accretech/examples/DF.mdf --dry-run --limit-rows 1
```

`run_wafer.py` is the outer prober loop: it parses an `.mdf` probe plan (`accretech_prober.parsers.read_mdf_probe_plan` → `(die_x, die_y)` per **PROB** die), and for each die calls `controller.move_to_die(x, y)` then the **same** `run_device()` `run_all.py` uses — with `run_info` derived from `settings/run_info.json` (position = die x/y, auto-incrementing `device_no`). The prober is the standalone `accretech-prober` package (under `accretech/`, wired as a **uv editable path dependency**; not modified by this repo). `--dry-run` swaps in a `NullProberController` so the whole loop runs with no GPIB. On error it calls `abort_and_safe_state()`; it always `separate()`s and closes the prober in `finally`.

Dispatch logic (how the runner is chosen)

Both `sweep.py` and `run_all.py` load the JSON and pick an engine based on which keys are present. **No engine is named in the config; it is inferred:**

Key present in JSON	Engine class	Spec loader
"lcr"	LCRRunner	load_lcr_spec
"series_resistance"	SeriesResistanceRunner	load_measurement_spec

Key present in JSON	Engine class	Spec loader
neither (has "smus")	SweepRunner	load_measurement_spec

All runners expose the **same interface**:

```
Runner(spec, dry_run, limit_rows, output_path).run() -> Path
```

That uniform contract is what lets `sweep.py/run_all.py` treat them interchangeably.

6. Why there are several engine files

The engines live in `core/` and exist because the three measurement families differ in **what happens inside the per-row loop** — the source/measure step — even though the surrounding scaffolding (load rows, connect, route matrix, shutdown) is identical. Rather than one giant branching loop, the variation is expressed through a small class hierarchy.

```

SweepRunner          (core/engine.py)          ← generic V/I sweep
└─ SeriesResistanceRunner (core/series_resistance_engine.py) ← overrides
   _run_row only

LCRRunner            (core/lcr_engine.py)      ← standalone, drives the 590

```

6.1 SweepRunner — `core/engine.py`

The generic engine and the base class. For each row it:

1. Routes the matrix to the row's crosspoints.
2. Walks the **Cartesian product** of all channels' sweep steps (`_build_combos`, built on `itertools.product` so the first channel is the outermost loop). At each combination it applies every channel's setpoint + hold in declaration order.
3. Measures **all** channels and appends **one CSV row per combination**, with column names generated from the channel roles (`Cathode Target V`, `Cathode Current A`, ...). Each row carries a flattened global `Step Index` plus a per-channel `<Label> Step Index`.

This covers guard leakage and dark current: source a voltage on one channel, hold the others, read the current. With only one multi-step channel the product collapses to that channel's steps, and the global `Step Index` equals its step index — which is what `report.py` relies on (see §6.4 / the note in `report.py`).

6.2 SeriesResistanceRunner — `core/series_resistance_engine.py`

Subclasses `SweepRunner` and **overrides only `_run_row`**. It reuses all the connection, matrix, and shutdown logic from the base class. The difference:

- It sources **two current** setpoints (the primary is in `"Current in A"` mode).
- It measures V & I at each, then computes $R_s = ((V_2 - V_1) + 2 \cdot \ln(I_2/I_1) \cdot 26\text{mV}) / (I_2 - I_1)$.

- It emits **one CSV row per pin** (the computed Rs), not one per step.

This is the textbook case for inheritance: same loop skeleton, different math at the measurement point.

6.3 LCRRunner — `core/lcr_engine.py`

A **standalone** runner (does *not* extend `SweepRunner`) because it drives an entirely different instrument (the 590, not the 4200) with a different port termination, a different configure/measure protocol, and no notion of multiple SMU channels. It keeps the same outer shape — load rows, connect, per-row loop, safe shutdown — and the same `.run()` -> `Path` contract, but the body sets a DC bias, triggers one C/G read, and records capacitance, conductance, and DC bias per pin.

It is a sibling rather than a subclass because almost nothing in the SMU-specific base would be reused; sharing it would mean fighting the inheritance rather than benefiting from it.

6.4 How the engine is "driven"

The engine is **driven by the spec dataclass**, which is in turn built from the JSON. The flow:

1. `sweep.py` reads the JSON, inspects its keys, and calls the matching `load*_spec()`.
2. The loader (`core/channel.py` or `core/lcr_channel.py`) parses the JSON into a **frozen dataclass** (`MeasurementSpec` / `LCRMeasurementSpec`). This is where parsing and any validation live. Nothing mutable or hardware-related is in the data model.
3. The chosen Runner is constructed with that spec and `.run()` is called. The runner reads everything it needs (addresses, channels, sweep profile) off the spec — it never re-reads the JSON.

So: **JSON** → **loader** → **immutable spec** → **runner**. The runner's behavior is fully determined by the spec it is handed; swapping configs swaps behavior without touching engine code.

7. Supporting modules

Module	Responsibility
<code>core/channel.py</code>	<code>MeasurementSpec/ChannelSpec/HardwareConfig/SweepStep</code> dataclasses + <code>load_measurement_spec()</code> (SMU configs). Channels keep JSON declaration order, which the engine uses as sweep-nest order.
<code>core/lcr_channel.py</code>	<code>LCRMeasurementSpec/LCRSpec</code> dataclasses + <code>load_lcr_spec()</code> (LCR config).
<code>core/port.py</code>	<code>PortWrapper</code> adapts a PyVISA resource to the <code>.write/.read/.query</code> API the SweepMe! drivers expect. <code>DryRunPort/DryRunResourceManager</code> stub hardware for <code>--dry-run</code> and print every command.
<code>switching_matrix.py</code>	<code>SwitchingMatrix707A</code> wraps the 707A driver (<code>apply_route</code> , <code>open_all</code> , <code>shutdown</code>); <code>normalize_matrix_config()</code> canonicalizes crosspoint strings.
<code>parameter_matrix.py</code>	Reads/normalizes the Excel workbook into row dicts.
<code>core/report.py</code>	<code>ReportWriter</code> turns the four measurement CSVs + <code>run_info.json</code> into the station's tab-separated <code>.txt</code> result with per-pin pass/fail bins and

Module	Responsibility
	limits.
<code>core/post_process.py</code>	SweepMe!-style post-processing script (standalone; not called by <code>sweep.py/run_all.py</code>).
<code>keithley4200.py</code> , <code>keithley590_LCR.py</code> , <code>keithley707a.py</code>	The actual SweepMe! instrument drivers (GPIB command protocols).
<code>helper/*.py</code>	Standalone GPIB diagnostics (<code>ping_4200.py</code> , <code>deep_flush.py</code> , <code>bus_reset.py</code> , <code>test_4200_*.py</code>). Run directly.

8. Outputs

- **sweep.py (single measurement)** — `results/<TestName>_<timestamp>.csv`, one row per sweep-nest combination (or per pin for series resistance / capacitance), columns derived from channel roles, plus a global `Step Index` and per-channel `<Label> Step Index`.
- **run_all.py / run_wafer.py (per device)** — a folder `results/wafer<id>/dev<NNN>_x<x>_y<y>/` holding the four measurement CSVs (`guard_leakage.csv`, `dark_current.csv`, `capacitance.csv`, `series_resistance.csv`) and the device report `result_dev<NNN>.txt`. `run_wafer.py` creates one such folder per probed die.
- **Report** — tab-separated `result_dev<NNN>.txt`, with a metadata header (from `run_info.json`), one line per pin-test with applied limits and a pass/fail bin, and a footer with total time, soft bin, and overall pass flag.

9. Environment

- Requires **Python 3.9.x exactly** (`requires-python = "==3.9.*"`).
- Install with `uv sync`.
- Instruments communicate over **GPIB via PyVISA**; the drivers are SweepMe! device classes that depend on `pysweepme`.
- Use `--dry-run` to exercise the full code path with no hardware attached — every instrument command is printed instead of sent.